




V300

LED Curing Light



5W Blue Light LED



Constant light output intensity

Not affected by battery power, ensuring stable curing performance.



Wavelength




400nm~515nm

Maximum Light Intensity




1200 mW/cm²


Mode 01:
Full Power Mode



Mode 02:
Progressing Mode




Mode 03:
Pulse Mode



Time Settings:

3s, 5s, 10s, 15s, 20s, 25s, 30s, 35s, 40s




SevereNote: Treatment timing should be determined depending on the severity of the patient's caries.



No disassemble


The built-in battery is capable of lasting over 400 full charge cycles.



Low battery alarm

The low battery alarm is indicated by a flashing display and an audible alert.

Ergonomic design enhances grip texture and comfort



M: Mode Selection

T: Time Setting

P: Power Button

01: Full Power Mode

02: Progressing Mode

03: Pulse Mode

